


# SHOEBOX for Hearing Services

A Mobile Solution for Out-of-the-Booth Testing, Wider Customer Reach & Lead Generation

*Advanced diagnostic capabilities, state-of-the-art lead generation, and cloud-based data management packed into one portable, easy-to-use solution.*




SHOEBOX® for Hearing Services is for the clinic owner ready to grow their business. We know it can be a challenge to attract new, qualified patients and prioritize appointment booking based on their hearing level. And we know that you probably wish to offer additional services outside of your clinic. If so, we can help. Both with increasing foot traffic and out-of-the-booth testing.

## SHOEBOX for Hearing Services:

- ✓ Self-administered screening
- ✓ Pure tone: air & bone with masking
- ✓ SRT/SDT testing
- ✓ Embedded inventory, surveys, or custom questionnaires
- ✓ Threshold tests automatically sent to Noah® 4
- ✓ Secure reporting and export using a HIPAA-compliant cloud portal

We offer a bundled solution of SHOEBOX QuickTest and SHOEBOX Audiometry Pro Edition. This highly portable, powerful combination of products offers testing functions on a single platform. Those interested in online hearing screening may also choose to add SHOEBOX Online.

# SHOEBOX for Hearing Services




SHOEBOX QuickTest can be set up in both attended and unattended kiosk mode at health fairs, in retail locations, waiting rooms, optometry offices, or almost any setting where your potential clients may gather. This quick, self-administered hearing screening check enables you to capture qualified leads in a fun and informative way.



When you need to perform threshold seeking tests, SHOEBOX Audiometry Pro offers a full range of audiometric capabilities including air and bone conduction with masking, speech reception threshold testing (SRT), speech discrimination testing, extended high frequency (EHF) testing, and sound generation ranging from -10dB to 120dB.



SHOEBOX Online enables participants to assess their hearing from any device - mobile phone, tablet, or laptop - and with the headphones they already own. SHOEBOX Online can help you reach and begin to engage a broader range of customers who are ready to take the next step in their hearing health journey.



The powerful combination of SHOEBOX QuickTest and SHOEBOX Audiometry offers a complete solution for hearing screening, lead generation, diagnostic evaluations, and tele-audiology testing to the mobile clinician. Those interested in online hearing screening may also choose to add SHOEBOX Online.