

























SHOEBOX Audiometry Symbols Legend: 4.x

Symbols that were in use prior to SHOEBOX 5.0










Audiogram Threshold Symbols

Symbol	Meaning	Symbol	Meaning	Symbol	Meaning
	Left air conduction		Right air conduction		Air
	No response to upper limits air conduction, left		No response to upper limits air conduction, right		Masked air
	Left masked air conduction		Right masked air conduction		Bone
	No response to upper limits masked air conduction, left		No response to upper limits masked air conduction, right		Too Noisy
	Unmasked bone, left side		Unmasked bone, right side		Unreliable
	No response to upper limits unmasked bone, left side		No response to upper limits unmasked bone, right side		Needs Masking
	Masked bone conduction left side		Masked bone conduction right side		Threshold should be re-tested
	No response masked bone left side		No response masked bone right side		No Response

Tabular View Abbreviations

Symbol	Meaning	Symbol	Meaning	Symbol	Meaning
NR	No response	M	Threshold was obtained using masking	TN	Result too noisy
+B	Bone testing suggested	+M	Masking suggested	?	Unreliable

Other Software Symbols

Symbol	Meaning	Symbol	Meaning	Symbol	Meaning
	Heard a sound		Did not hear a sound		New Patient
	Patient Information		Upload to web portal		Notes
	Add Item / Create New		Reorder Item		Warning