


SHOEBOX QuickTest

A Simple Solution for New Customer Acquisition

An easy-to-use hearing screening test to help you find qualified new customers for your hearing clinic. Qualified leads for your hearing services clinic, at a fraction of the cost of traditional marketing.



- ✓ Configurable testing options - to suit your specific needs
- ✓ Custom branding - to promote your company and maintain consistency
- ✓ Supports 19 languages, HIPPA Compliant
- ✓ Ability to manage a fleet of iPads through a web portal

We know it's a challenge to attract new, qualified patients and prioritize appointments based on their hearing levels. We believe you deserve an easy-to-use and cost-effective solution. SHOEBOS® QuickTest can be used in both attended and unattended kiosk mode in many settings where your potential clients may gather. SHOEBOS QuickTest can help you identify your most qualified patients, sort them into different hearing levels, and recommend the next step in their hearing health journey.

With automatic syncing of results to our cloud-based portal, you can quickly and easily identify your best candidates for scheduling follow-up appointments—no additional work required.

SHOEBOS QuickTest is the ultimate companion to our SHOEBOS Pro edition, which can be used to perform diagnostic hearing testing after using SHOEBOS QuickTest.




Easy-to-use

Simple, self-administered hearing screening test



Qualified Leads

Prioritize follow-up based on participants' results



Quick Triage

Guide participants to the right service based on their hearing levels

SHOEBOX QuickTest Features

Equipment Compatible with most iPads, see full list

Supported Transducers RadioEar DD450

Functionality

- Configure up to 5 frequencies & choose between 2 or 3 hearing levels
- Create custom action pages to inform and guide participants on your recommended next steps
- Customize to your brand images, logo and color
- Available for participants in 19 different languages
- Enable multi-language feature for participants to select from
- Give participants the choice to email themselves their results and follow up with you directly
- Optional request a booking feature to prioritize those ready to book today, admin notification feature for a daily summary
- Turn on our questionnaire to capture more insights
- Admin feature to review results together in assisted mode
- Participant overview to edit and review post test
- Collect contact information and tailor to your specific requirements
- Set your own privacy policy to conform to local laws and regulations
- Centralized command and control for all devices
- Manage and analyze screening results and additional information from all devices via a web portal

Sound Generation Range: 500Hz - 8000Hz
0 dB(HL) to 90 dB(HL)

Subscription Includes:

- ✓ Software and license to use on one iPad device
- ✓ Cloud-based web portal
- ✓ One hour training/setup session

Additional Options:

- ✓ Kiosk floor stand or desktop stand bundle
- ✓ Remote iPad Management

Hardware Requirements:

- ✓ Compatible iPad
- ✓ Compatible Transducers